

Product Name :
Recovery Bed**Product Code :**
NMCBS106006**Description :**

Recovery Bed We are among the reputed organizations, engaged in offering our clients with Recovery Bed. All the beds are designed in compliance with international quality norms & standards, using premium raw material, which is sourced from the reliable and certified vendors of the market. In addition to this, we have a team of expert quality controllers, which thoroughly checks all the products on various parameters to ensure flawlessness. Other Information: $\frac{1}{2}$ Overall Size: L 2010 mm x W 865 mm x H 590 mm $\frac{1}{2}$ Frame: 50 mm x 25 mm $\frac{1}{2}$ CRCA Rectangular tube $\frac{1}{2}$ Two section perforated CRCA sheet top $\frac{1}{2}$ 2 Separate screws for back rest, knee rest trendelenburg & reverse trendelenburg position by attached folding handles $\frac{1}{2}$ Detachable sun-mica paneled bows $\frac{1}{2}$ 125 mm dia heavy duty wheels; 2 with brakes $\frac{1}{2}$ 4 flush I.V. rod locations $\frac{1}{2}$ Pre-treated and epoxy powder coated Recovery Bed, Hospital Bed exporters, Hospital Bed manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: [SITE](#)