

Product Name :
Recovery Bed

Product Code :
NMCBS106006



Description :

Recovery Bed We are among the reputed organizations, engaged in offering our clients with Recovery Bed. All the beds are designed in compliance with international quality norms & standards, using premium raw material, which is sourced from the reliable and certified vendors of the market. In addition to this, we have a team of expert quality controllers, which thoroughly checks all the products on various parameters to ensure flawlessness. Other Information: • Overall Size: L 2010 mm x W 865 mm x H 590 mm • Frame: 50 mm x 25 mm • CRCA Rectangular tube • Two section perforated CRCA sheet top • 2 Separate screws for back rest, knee rest trendelenburg & reverse trendelenburg position by attached folding handles • Detachable sun-mica paneled bows • 125 mm dia heavy duty wheels; 2 with brakes • 4 flush I.V. rod locations • Pre-treated and epoxy powder coated Recovery Bed, Hospital Bed exporters, Hospital Bed manufacturers

NAUGRA[®]

Medical Hospital Equipments Manufacturers

Website:Medical Hospital Equipments Manufacturers,