

Product Name :
Medical Fowler Bed Manual**Product Code :**
NMCBS106014**Description :**

Medical Fowler Bed Manual

Technical Specification :

Medical Fowler Bed Manual with mechanically operated back rest and knee rest by super-smooth crank mechanism.

Size: 2140 940 540 mm Removable & interchangeable high quality beautiful ABS engineering plastic head panel and foot panel. Collapsible Stainless Steel side railings. Epoxy coated mild steel frame work and 4 section perforated top. 125mm dia noiseless castors (all with individual brakes). Provision for I.V. Rod. Elit Fowler Bed Manual, Hospital Bed exporters, Hospital Bed manufacturers