

**Product Name :**  
Orthopaedic Bed Two Sectioned Adjustable Knee Rest**Product Code :**  
NMCBS106017**Description :**

Orthopaedic Bed Two Sectioned Adjustable Knee Rest is a Orthopaedic bed with mechanically operated back rest & knee rest all by crank mechanism. Knee-rest portion sectioned in two parts for patient's convenience. Size: 2030 x 900 x 500 mm. Epoxy coated mild steel 5 section top. Frame work made of stainless steel. Balkan Frame with monkey bars and monkey chains. Pulley attachment (without weight chain and weight hangers). Mounted on protective stumps. Provision for I.V. Rod. Freight saving knock-down construction. Type:- Without castor With castor Orthopaedic Bed Two Sectioned Adjustable Knee Rest, Hospital Bed exporters, Hospital Bed manufacturers