



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Hospital Fowler Beds Mechanically

Product Code :
NMCBS106019



Description :

Hospital Fowler Beds Mechanically Fowler Beds Mechanically Fowler Bed with various body positions. Fowler Beds Mechanically Backrest & Knee rest section maneuvered by two separate screws Fowler Beds Mechanically Framework strong & precise rectangular steel tubes mounted on 10 cm dia Fowler Beds Mechanically Swiveling Castors two with brake. Fowler Beds Mechanically Tubular head and foot bows unequal heights with vertical tubes supports. Fowler Beds Mechanically Provisions for four saline rod locations Fowler Beds Mechanically Overall Size: 198L x 90W x 60 H (Approx) cms Finish: Pre-treated to prevent rust and Epoxy Powder Coated. Hospital Fowler Beds Mechanically, Hospital Bed exporters, Hospital Bed manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,