

Product Name :
Electrical Intensive Care Bed**Product Code :**
NMCBS106022**Description :**

Electrical Intensive Care Bed Intensive care Hospital Bed for comfortable recovery of patient under cumbersome conditions. ½ Four section CRCA sheet Top, uniformly perforated fixed to the frame work & base made of ERW Steel tube /rectangular tubes ½ MS channels mounted on four special Swivelling castors with 2 brakes ½ full length collapsible stainless steel safety railing. ½ Adjustable telescopic I.V.Rod. Remote control system offers:

½ : Electrical Intensive Care Bed Is High & low position ½ : Electrical Intensive Care Bed is Headrest position ½ : Electrical Intensive Care Bed Trend lenburg/reverse trend lenburg position ½ : Electrical Intensive Care Bed Size: 198(L) x 90(W) x 60 to 80(H) cms ½ Finish: Pre-treated to prevent rust and Epoxy Powder Coated.

Electrical Intensive Care Bed, Hospital Bed exporters, Hospital Bed manufacturers