

Product Name :
Electrical Intensive Care Bed (With Remote)

Product Code :
NMCBS106023



Description :

Electrical Intensive Care Bed (With Remote) Electrical Intensive Care Bed (With remote) Intensive care Hospital Bed for comfortable recovery of patient under cumbersome conditions. Electrical Intensive Care Bed (With remote) Four section CRCA sheet Top, uniformly perforated fixed to the frame work & base made of ERW Steel tube /rectangular tubes Electrical Intensive Care Bed (With remote) MS channels mounted on four special Swivelling castors with 2 brakes Electrical Intensive Care Bed (With remote) full length collapsible stainless steel safety railing. Electrical Intensive Care Bed (With remote) Adjustable telescopic I.V.Rod. Remote control system offers: Electrical Intensive Care Bed (With remote) High & low position Electrical Intensive Care Bed (With remote) Headrest position Electrical Intensive Care Bed (With remote) Trendlenburg/reverse trendlenburg position Electrical Intensive Care Bed (With remote) Size: 198(L) x 90(W) x 60 to 80(H) cms Electrical Intensive Care Bed (With remote) Finish: Pre-treated to prevent rust and Epoxy Powder Coated. . Electrical Intensive Care Bed (With Remote), Hospital Bed exporters, Hospital Bed manufacturers