



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
Monitor Trolley For Scopic Surgery

Product Code :
NMCBS109013



Description :

Monitor Trolley For Scopic Surgery Overall size W 660 mm x D 560 mm x H 1498 mm Features: Available in full Stainless steel & epoxy powder coated Adjustable Shelves 125 mm dia. heavy duty casters with brakes Monitor Trolley For Scopic Surgery, Hospital Trolley exporters, Hospital Trolley manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE