

**Product Name :**  
Mayo Type Trolley Height Adjustable**Product Code :**  
NMCBS109017**Description :**

Mayo Type Trolley Height Adjustable is a Trolley with Stainless Steel framework. It is Used in the O.T. rooms for keeping instruments during surgery. Size: 600 x 400 x 850 / 1000 mm. Adjustable height from 850 - 1000 mm. 2 Swivel castors (with brakes). Mayo Type Trolley Height Adjustable, Hospital Trolley exporters, Hospital Trolley manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [SITE](http://SITE)