

**Product Name :**  
Hydraulic Emergency Trolley**Product Code :**  
NMCBS109020**Description :**

Hydraulic Emergency Trolley Hydraulic Emergency Trolley \* Two section top. \* Height adjusted by foot operated imported hydraulic pump from 660 mm to 910 mm. \* Gas spring assisted Trendelenburg/Reverse Trendelenburg positions \* Rubber covered handles, accessories tray, oxygen cylinder holder, SS telescopic IV rod and swing away SS side rails \* Pre-treated and powder coated. \* Four swivel castors, 20cms dia, two with total lock. \* X-ray Permeable removable stretcher, backrest raised on ratchet. \* Overall approx size: 185L x 60W x 65-95H (CM)  
Hydraulic Emergency Trolley, Hospital Trolley exporters, Hospital Trolley manufacturers

**NAUGRA**<sup>®</sup>**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,