

Product Name :
Hydraulic Emergency Trolley**Product Code :**
NMCBS109020**Description :**

Hydraulic Emergency Trolley Hydraulic Emergency Trolley * Two section top. * Height adjusted by foot operated imported hydraulic pump from 660 mm to 910 mm. * Gas spring assisted Trendelenburg/Reverse Trendelenburg positions * Rubber covered handles, accessories tray, oxygen cylinder holder, SS telescopic IV rod and swing away SS side rails * Pre-treated and powder coated. * Four swivel castors, 20cms dia, two with total lock. * X-ray Permeable removable stretcher, backrest raised on ratchet. * Overall approx size: 185L x 60W x 65-95H (CM)

Hydraulic Emergency Trolley, Hospital Trolley exporters, Hospital Trolley manufacturers