

Product Name :
Emergency & Recovery Trolleys**Product Code :**
NMCBS109027**Description :**

Emergency & Recovery Trolleys
Emergency & Recovery Trolley Made of heavy duty steel tubes.
Removable X-Ray translucent provided with a tray for holding X-ray cassettes underneath the entire length of the stretcher.
Emergency & Recovery Trolley Height adjustment by hydraulic system, trendelenburg and reverse trendelenburg by pressure cylinder.
Emergency & Recovery Trolley Fitted with 15cms dia castors- two with brakes and rubber buffers on all sides.
Emergency & Recovery Trolley Provided with SS swing away railing, SS telescopic saline rod, storage tray and cylinder cage.
Emergency & Recovery Trolley Size: 210L x 65W x 65 to 95 H cms.
Finish: Pre-treated to prevent rust & Epoxy powder coated.
Emergency & Recovery Trolleys, Hospital Trolley exporters, Hospital Trolley manufacturers