

**Product Name :**  
Hospital Locker**Product Code :**  
NMCBS110003**Description :**

Hospital Lockers Backed by our efficient workforce, we are able to manufacture and supply a wide range of Hospital Lockers. Our range encompasses Multipurpose Bed Side Lockers, Bed Side Locker (Deluxe), Bed Side Locker, Instrument Cabinet, Bedside Cabinet ABS and Bedside Cabinet ABS. All our products are designed and manufactured using premium raw material, which is sourced from the reliable vendors of the market. In addition to this, we have the capability to customize all the products as per the specifications laid down by the customers. Specifications:  $\frac{1}{2}$  Overall Size: W 410mm X D 370mm X H 900mm  $\frac{1}{2}$  Size after packing: W 840mm X D 500mm X H 190mm Hospital Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**Website: [MedicalHospitalEquipmentsManufacturers.com](http://MedicalHospitalEquipmentsManufacturers.com)