



Email : sales@naugramedical.com
Phone: 9896600003

Product Name :
ABS Bed Side Lockers

Product Code :
NMCBS110008



Description :

ABS Bed Side Lockers The Hospital ABS Bedside Locker is equipped with a locker machine, pressed from CRCA sheet, with all three sides closed, having one box and a drawer. The Hospital Bedside Locker has a stainless steel top and mounted on two swivel castors. ABS Bed Side Lockers, Hospital Lockers exporters, Hospital Lockers manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE