

**Product Name :**  
Ceiling Lamp with 4 Reflector**Product Code :**  
NMCBS113002**Description :**

Ceiling Lamp with 4 Reflector & Dome Head: 515mm Dia. & Halogen Bulbs :4nos. 12v x 50 watt & Light intensity at 1 mt.: 70,000 Lux Intensity Control: 4 Stage & Height Adjustment: 550 mm & Action Radius: 1550mm & Possible Movements: Radial Angular & Axial Color & Temperature: 4200 K & 300K & Temp.rise in field: 3-6 c from Amb.Temp & Low voltage unit with CVT: 1 No. (C.V.T optional) & Input Supply Voltage: 220 Volt A.C & 10% 50Hz & 2% Ceiling Lamp with 4 Reflector, Operation Theater Lights exporters, Operation Theater Lights manufacturers

**NAUGRA®****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,