

**Product Name :**  
Operation White Light 7 Reflector**Product Code :**  
NMCBS113004**Description :**

Operation White Light 7 Reflector We have gained expertise in providing our customers with Operation White Light 7 Reflector. These products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and excellent performance. In addition to this, we pack all the products using premium packaging material to ensure complete safety during transportation. Dimensions: Dome Head:715mm Dia. Halogen Bulbs:7nos. 12v x 50 watt Light intensity at 1 mt.:1,20,000 Lux Intensity Control:4 Stage Height Adjustment:500mm Action Radius:1850mm Possible Movements:Radial Angular & Axial Color Temperature:4200 K  $\pm$  300K Temp rise in field:3 $\frac{1}{2}$  - 6 $\frac{1}{2}$  c from Amb. Temp Low voltage unit with CVT:1 No. (C.V.T optional) Input Supply Voltage:220voltA.C  $\pm$  10% 50Hz  $\pm$  2% Operation White Light 7 Reflector, Operation Theater Lights exporters, Operation Theater Lights manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website:Medical Hospital Equipments Manufacturers,