

Product Name :
Anti Bedsore Mattress**Product Code :**
NMCBS115002**Description :**

Anti Bedsore Mattress ĩ½ Anti Bedsore Mattress based on alternative pressure pad & pump system ĩ½ Power supply: 220V/50Hz, 1A/250V, Air output: 6-8 liter/min ĩ½ Pressure range: 30 mmHg ĩ½ 100 mmHg ĩ½ Top cover: Nylon / PU with button ĩ½ Size:190cm(L) x 85cm(W) x 11cm (H) ĩ½ Exchangeable air cell, water proof. Low Noisy, Long life time benefit Size: ĩ½ Without CPR valve ĩ½ With CPR valve Anti Bedsore Mattress, Hospital Beds Mattress exporters, Hospital Beds Mattress manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**Website: SITE