

Product Name :
Anti Bedsore Mattress**Product Code :**
NMCBS115002**Description :**

Anti Bedsore Mattress i.e. Anti Bedsore Mattress based on alternative pressure pad & pump system i.e. Power supply: 220V/50Hz, 1A/250V, Air output: 6-8 liter/min i.e. Pressure range: 30 mmHg i.e. 100 mmHg i.e. Top cover: Nylon / PU with button i.e. Size: 190cm(L) x 85cm(W) x 11cm (H) i.e. Exchangeable air cell, water proof. Low Noisy, Long life time benefit Size: i.e. Without CPR valve i.e. With CPR valve Anti Bedsore Mattress, Hospital Beds Mattress exporters, Hospital Beds Mattress manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,