



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Bedside Screen Superior

Product Code :
NMCBS117002



Description :

Bedside Screen Superior We have established ourselves as a renowned manufacturer and supplier of Bedside Screen Superior. All our products are widely acclaimed by the customers for their corrosion & abrasion resistance, sturdy construction, low maintenance and easy handling. These products are designed and developed using premium raw material, which is sourced from the reliable vendors of the market. Specifications: $\frac{1}{2}$ Size: 168 H x 244W cms $\frac{1}{2}$ Supplied with Rod & Curtains $\frac{1}{2}$ Special Knockdown construction $\frac{1}{2}$ Frame work made of Epoxy Coated Mild Steel with castors Type: $\frac{1}{2}$ Bedside Screen - Stainless Steel - 3 folds $\frac{1}{2}$ Bedside Screen - Epoxy Coated Mild Steel 4 folds $\frac{1}{2}$ Bedside Screen - Epoxy Coated Mild Steel - 3 folds $\frac{1}{2}$ Bedside Screen - Stainless Steel - 4 folds Bedside Screen Superior, Hospital Bedside Screens exporters, Hospital Bedside Screens manufacturers



Medical Hospital Equipments Manufacturers

Website: Medical Hospital Equipments Manufacturers,