

Product Name :
Labor & Delivery Bed with Accessories**Product Code :**
NMCBS120005**Description :**

Labor & Delivery Bed with Accessories Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes. Leg section can be completely recessed under the body section
Size :
Body section: 1000(L) x 800(W) x 750(H)mm., Leg section: 900(L) x 800(W) x 750(H)mm
Supports for leg holders with solid mechanism to fix leg holders at required height and width
Mattress-pad 2 sections
Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width
Epoxy coated tubular steel frame
Mattress-pad density : 21-25 kg/m³
Carrying capacity: 150 - 180 kg
Mattress-pad: Approx. 60(H)mm
Supplied with S.S. bowl
Knockdown construction Labor & Delivery Bed with Accessories, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers

NAUGRA®**Medical Hospital Equipments Manufacturers**Website: SITE