

Product Name :
Dismantled Gynecology Table**Product Code :**
NMCBS120006**Description :**

Dismantled Gynecology Table Made in 3 sections, adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips
Size: Table 1800L x 600W x 800H mm. Tubes : 30 x 1.5 mm
Supports for leg holders with solid mechanism to fix leg holders at required height and width
Support ring at proximity of seat section to support round basin-bowl
Supplied with 2 leg holders with canvas straps adjustable in height and width
Epoxy coated tubular steel frame
Tubes: 30 x 1.5 mm
Mattress-pad: density: 21 - 25 kg/m³
Carrying capacity: 150-180 kg
Mattress-pad: 30H mm
Knockdown construction
Supplied with 4 liter S.S. Bowl Dismantled Gynecology Table, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers

NAUGRA[®]**Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,