

**Product Name :**  
Dismantled Gynecology Table**Product Code :**  
NMCBS120006**Description :**

Dismantled Gynecology Table Made in 3 sections, adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips  
Size: Table 1800L x 600W x 800H mm. Tubes : 30 x 1.5 mm  
Supports for leg holders with solid mechanism to fix leg holders at required height and width  
Support ring at proximity of seat section to support round basin-bowl  
Supplied with 2 leg holders with canvas straps adjustable in height and width  
Epoxy coated tubular steel frame  
Tubes: 30 x 1.5 mm  
Mattress-pad: density: 21 - 25 kg/m<sup>3</sup>  
Carrying capacity: 150-180 kg  
Mattress-pad: 30H mm  
Knockdown construction  
Supplied with 4 liter S.S. Bowl Dismantled Gynecology Table, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers