



Email : sales@naugramedical.com

Phone: 9896600003

Product Name :
Delivery Bed

Product Code :
NMCBS120007



Description :

Delivery Bed Size: 1800L x 550W x 800H mm Adjustable back rest and leg section with lever mechanism to varying positions, capable of tilting for obstetric and gynecological use. Legs mounted on rubber tips Supports for leg holders with solid mechanism to fix leg holders at required height and width Supplied with 2 knee crutches, adjustable in height and width Framework made of mild steel epoxy coated finish with upholstered top Supplied with S.S. bowl Knockdown construction Delivery Bed, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers



Medical Hospital Equipments Manufacturers

Website: SITE