

**Product Name :**  
Labor & Delivery Bed**Product Code :**  
NMCBS120008**Description :**

Labor & Delivery Bed Made in 2 sections, body section with legs mounted on rubber tips, legs section on swivel castors with brakes Size : Body section: 1000L x 800W x 750H mm., Leg section: 900L x 800W x 750H mm Supports for leg holders with solid mechanism to fix leg holders at required height and width, Supplied with 2 legs holders with canvas straps OR 2 knee crutches (as desired) adjustable in height and width. Epoxy coated tubular steel frame Carrying capacity: 150 - 180 kg Supplied with S.S. bowl Labor & Delivery Bed, Hospital Obstetric Tables exporters, Hospital Obstetric Tables manufacturers

**NAUGRA<sup>®</sup>****Medical Hospital Equipments Manufacturers**

Website: Medical Hospital Equipments Manufacturers,