

Product Name :
Material Storage Cabinet

Product Code :
NMCBS124007



Description :

Material Storage Cabinet Material Storage Cabinet. Owing to our rich industrial experience, we are offering Material Storage Cabinet that are used in clean rooms for storage of material. Widely demanded by various pharmaceutical, hospitals and O/T and various electronic industries, our products are featured with U.V. Germicidal light, HEPA & Pre-filter, motor blower, static pressure manometer, fluorescent light and shelves. In addition to this, we provide customization as per the specifications laid down by the customers. Features are given below: Material Storage Cabinet are Large space Material Storage Cabinet are Safe storage Material Storage Cabinet are Long durability Specifications: Working Size 3' 1/2 x 2' 1/2 x 5' 1/2 4' 1/2 x 2' 1/2 x 5' 1/2 Size of HEPA Filters 3' 1/2 x 1' 1/2 x 6' 1/2 4' 1/2 x 1' 1/2 x 6' 1/2 No. of HEPA Filters 1 1 No. of Pre Filter 1 2 No. of Pre Filter 3 3 Whole Cabinet made of commercial board & laminated with sun mica Whole Cabinet made of cold rolled steel sheet duly powder coated inside & outside for longer life durability Material Storage Cabinet, Medical Cabinets Cupboards exporters, Medical Cabinets Cupboards manufacturers